

InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody

Catalog No: #V1031

Package Size: #V1031-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1739
Isotype	IgG1, lambda
Purification	Protein A/G purified from cell culture supernatant.
Applications	Blocking, ELISA
Species Reactivity	Influenza A virus (A/Brisbane/59/2007(H1N1)) (D5F1Q8), Influenza A virus (A/Viet Nam/1203/2004(H5N1)) (Q6DQ33), Influenza A virus (strain A/Hong Kong/1/1968 H3N2) (Q91MA7), Influenza A virus (A/Netherlands/219/2003(H7N7)) (Q6VMK1)
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	B4URD6, A8C8J4, Q289M7
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.