

InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody

Catalog No: #V1033

Package Size: #V1033-1 100ug



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

|                    |                                                                                                                                                                           |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody                                                                                                                |
| Host Species       | Human                                                                                                                                                                     |
| Clonality          | Monoclonal                                                                                                                                                                |
| Clone No.          | VS2088                                                                                                                                                                    |
| Isotype            | IgG1, kappa                                                                                                                                                               |
| Purification       | Protein A/G purified from cell culture supernatant.                                                                                                                       |
| Applications       | ELISA, Neutralization                                                                                                                                                     |
| Species Reactivity | Influenza A virus (strain A/Hong Kong/212/2003 H5N1 genotype Z+)                                                                                                          |
| Conjugates         | Unconjugated                                                                                                                                                              |
| Target Name        | Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA                                                                                                       |
| Accession No.      | Q6J8F6                                                                                                                                                                    |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                          |
| Storage            | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.