

Anti-Influenza A virus (H3N2) HA/Hemagglutinin Recombinant Antibody

Catalog No: #V1036

Package Size: #V1036-1 50ug #V1036-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Influenza A virus (H3N2) HA/Hemagglutinin Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1030
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Influenza A virus (H3N2)
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	P03435
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.


Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000


Images



Detects recombinant H3N2 HA in indirect ELISA



SDS-PAGE for Anti- $\square$ Influenza A virus HA/Hemagglutinin Antibody



Western blot analysis was performed using anti- $\square$ monoclonal antibody at 1ug/mL on various samples. Lane 1: ATGE7523 200ng Lane 2: 1* loading buffer

Note: This product is for in vitro research use only and is not intended for use in humans or animals.