

Anti-Influenza A virus HA/Hemagglutinin Recombinant Antibody

Catalog No: #V1039



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1039-1 50ug #V1039-2 100ug


Description

Product Name	Anti-Influenza A virus HA/Hemagglutinin Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1614
Isotype	IgG1,kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Influenza A virus
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	B4URD6
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.


Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



Detects recombinant H1N1 HA in indirect ELISA



SDS-PAGE for Anti-Influenza A virus HA/Hemagglutinin Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.