

Anti-Influenza A virus (H5N1) HA/Hemagglutinin Recombinant Antibody

Catalog No: #V1042

Package Size: #V1042-1 50ug #V1042-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Influenza A virus (H5N1) HA/Hemagglutinin Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1667
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Influenza A virus (strain A/Vietnam/1203/2004 H5N1)
Conjugates	Unconjugated
Target Name	HA, Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain
Accession No.	A0A076G8X7
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.