

Anti-Influenza A virus NP/Nucleoprotein Recombinant Antibody

Catalog No: #V1043



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1043-1 50ug #V1043-2 100ug


Description

Product Name	Anti-Influenza A virus NP/Nucleoprotein Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1816
Isotype	VHH-mFc
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM, IF
Species Reactivity	Influenza A virus (strain A/Wilson-Smith/1933 H1N1) (Influenza A virus (strain A/WS/1933 H1N1)), Influenza A virus (A/WSN/1933(H1N1))
Conjugates	Unconjugated
Target Name	Nucleoprotein, Nucleocapsid protein, Protein N, NP
Accession No.	P15682
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A virus NP/Nucleoprotein Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.