

Anti-Influenza A virus RNA polymerase complex Recombinant Nanobody

Catalog No: #V1048

Package Size: #V1048-1 50ug #V1048-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-Influenza A virus RNA polymerase complex Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	VS1645
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA
Species Reactivity	Influenza A virus (strain A/Northern Territory/60/1968 H3N2) (Influenza A virus (strain NT60)) (Influenza A virus (strain A/NT/60/1968 H3N2)) (P03434 & P03432), Influenza A virus (A/Northern Territory/60/1968(H3N2)) (Q1AR05 & Q910D6)
Conjugates	Unconjugated
Target Name	Polymerase acidic protein, 3.1.-., RNA-directed RNA polymerase subunit P2, PA
Accession No.	P03434
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A virus RNA polymerase complex Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.