

Anti-Siltuximab Neutralizing Monoclonal Antibody

Catalog No: #V2742




Package Size: #V2742-1 50ug #V2742-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Siltuximab Neutralizing Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AD2048
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Immunogen Description	E. coli - derived recombinant protein.
Conjugates	Unconjugated
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Images



SDS-PAGE for Anti-Siltuximab Neutralizing Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.