

Anti-NeuGcGM3/GM3 Recombinant Antibody

Catalog No: #V8001



Package Size: #V8001-1 50ug #V8001-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-NeuGcGM3/GM3 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	GL2029
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IHC, WB
Species Reactivity	General
Conjugates	Unconjugated
Target Name	NeuGcGM3, Ganglioside N-glycolyl GM3, N-glycolylneuraminic acid -NeuGc, NGNA-GM3 ganglioside, NeuGcGM3
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SEC-HPLC detection for Anti-NeuGcGM3/GM3 Antibody .



SDS-PAGE for Anti-NeuGcGM3/GM3 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.