

Anti-TAG-72 Reference Antibody

Catalog No: #V8006

Package Size: #V8006-1 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-TAG-72 Reference Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2012
Isotype	IgG1-nd
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Bioactivity:B FACS, FunctionalB assay, ResearchB inB vivo
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	TAG-72
Accession No.	TAG-72
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Research Grade HuCC49.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.