

Anti-Ganglioside GM2 Recombinant Antibody

Catalog No: #V8027



Package Size: #V8027-1 50ug #V8027-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Ganglioside GM2 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	GL2039
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IHC
Species Reactivity	General
Conjugates	Unconjugated
Target Name	Ganglioside GM2, Monosialoganglioside GM2, GM2 lipid
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SEC-HPLC detection for Anti-Ganglioside GM2 Antibody .



SDS-PAGE for Anti-Ganglioside GM2 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.