

Anti-Human TAG72 Recombinant Antibody

Catalog No: #V8033



Package Size: #V8033-1 50ug #V8033-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Human TAG72 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	GL2001
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM, IP, RIA, WB
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	TAG-72, TAG72, HMW mucin-like glycoprotein, Cancer antigen 72-4, CA 72-4, CA 72-4/TAG72, STn, TAG-72
Accession No.	TAG-72
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Human CA72-4/TAG72 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.