

InVivoMAb Anti-sLea/CA19-9 Antibody

Catalog No: #V8035



Package Size: #V8035-1 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-sLea/CA19-9 Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2088
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	CA19-9, CA199, carbohydrate antigen sialyl Lewis(a), sLea, CA19.9, Sialyl Lewis A
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-sLea/CA19-9 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.