

Anti-sLea/CA19-9 Recombinant Antibody

Catalog No: #V8036



Package Size: #V8036-1 50ug #V8036-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-sLea/CA19-9 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2016
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	Blocking, ELISA, FCM, IF, IHC
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	CA19-9, CA199, Carbohydrate antigen199, Carbohydrate antigen 19-9, Sialyl Lewis A
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-sLea/CA19-9 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.