

Anti-Ganglioside GD3 Recombinant Antibody

Catalog No: #V8041



Package Size: #V8041-1 50ug #V8041-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Ganglioside GD3 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2032
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	General
Conjugates	Unconjugated
Target Name	Disialoganglioside GD3, GD3 ganglioside, α -N-acetylneuraminosyl-(2 β 8)- α -N-acetylneuraminosyl-(2 β 3)- β -D-galactosyl-(1 β 4)- β -D-glucosyl-(1 β 1')-ceramide
Accession No.	ChEBI: 28424, PubChem CID: 20057323
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Ganglioside GD3 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.