

InVivoMAb Anti-Globo H Antibody

Catalog No: #V8051



Package Size: #V8051-1 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-Globo H Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	GL2080
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FuncS
Conjugates	Unconjugated
Target Name	Blood group H antigen
Accession No.	ChEBI: 62361
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-Globo?H Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.