

InVivoMAb Anti-Human SSEA-4 Antibody

Catalog No: #V8053



Package Size: #V8053-1 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-Human SSEA-4 Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	GL2082
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FuncS
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	SSEA4, OB1-898, OB1 898, OB1898, OBI-898, OBI 898, OBI898
Accession No.	ChEBI: 141011
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-SSEA4 antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.