

Anti-Lewis Y antigen Recombinant Antibody

Catalog No: #V8057



Package Size: #V8057-1 50ug #V8057-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Lewis Y antigen Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2076
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	Alpha-(1, 3)-fucosyltransferase 4, 4-galactosyl-N-acetylglucosaminide 3-alpha-L-fucosyltransferase, 2.4.1.152, ELAM-1 ligand fucosyltransferase, Fucosyltransferase 4, Fucosyltransferase IV, Fuc-TIV, FucT-IV, Galactoside 3-L-fucosyltransferase, FUT4, ELFT, FCT3A
Accession No.	ChEBI: 62287
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Lewis Y antigen Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.