

Anti-Human CD15/SSEA-1/Lewis x Recombinant Antibody , PerCP

Catalog No: #V8074

Package Size: #V8074-1 50T #V8074-2 100T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-Human CD15/SSEA-1/Lewis x Recombinant Antibody , PerCP
Host Species	Mouse
Clonality	Monoclonal
Clone No.	GL2092
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM
Species Reactivity	Human
Conjugates	PerCP
Target Name	Alpha-(1, 3)-fucosyltransferase 4, ELAM-1 ligand fucosyltransferase, Fucosyltransferase 4, Fucosyltransferase IV, Galactoside 3-L-fucosyltransferase, Fuc-TIV, FCT3A, FUT4, ELFT, FucT-IV
Accession No.	P22083
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.

Application Details

Purity:>90% as determined by SDS-PAGE.

Images



Flow-cytometry using PerCP anti-human CD15 antibody. THP-1 cells were fixed and permeabilized, then stained with an irrelevant antibody (Blue Histogram) or an PerCP anti-human CD15 monoclonal antibody (Yellow Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, and cells analysed on a NovoCyt Flow Cytometer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.